

Styletto.

The world's first rechargeable SLIM RIC that combines style with high-tech hearing.



Styletto Features & Benefits

Exceptional style

SLIM RIC hearing instrument you wear with confidence

Portable rechargeability


Recharging case for easy handling and long wearing time

A natural hearing experience

Signia Nx for the most natural own voice and clear sound in every situation

Practice differentiation

Offer the proven solution that patients prefer and remember



touchControl App and miniPocket

To enhance its discretion, Styletto is remotely controllable via the touchControl App and the miniPocket.

- Contemporary design with clean surfaces
- Portable rechargeability
 - 3 hour full charge resulting in up to 19 hours daily use
 - 30 minute fast charge resulting in 5 additional hours use
 - 3 full charges on the go
- Most natural own voice and best speech understanding in noise with Signia Nx platform
- Robust fixed receiver and IP68 rated[†] for unmatched reliability
- Size M receiver (60/119) available in lengths 2 & 3



A new service concept that provides convenience for you and your patients. Request a replacement Styletto to be received prior to sending in original device.

Performance levels

7

5


3

sDemo

Uncompromising Audiological Performance - Signia Nx technology is **clinically-proven** to improve the most common communication issues.

- 80% improvement in acceptance of a wearer's own voice with Own Voice Processing*
- 25% better speech recognition in noisy environments for Signia wearers over people with normal hearing with Narrow-Directionality (offer your patient "better than normal hearing**")
- 60% reduced listening effort in demanding environments lessening wearer fatigue*

Type	Styletto
Performance levels	7 / 5 / 3



Product details	
Matrix (2 ccm)	M 60/119 dB
Battery size	Li-Ion power cell
Power consumption	1.4 mA

Features	
Frequency channels/ Compression controls	7Nx - 48/20 5Nx - 32/16 3Nx - 24/12
Hearing programs	6
Ultra HD e2e wireless	●
Direct streaming	-
TwinPhone	● [†]
TeleCare 2.0	●
Tinnitus therapies	●
CROS/BiCROS	-
Push button	-
Volume control	-
Rechargeable	●
T-coil	-
Battery door on/off function	-
Audio input	-

Modifications	
Color conversion kit	-

App/Accessory compatibility	
myControl App	-
myHearing App	●
touchControl App	●
StreamLine TV	-
StreamLine Mic	-
miniPocket	●
Noahlink™ Wireless	-
ConnexxAir/ConnexxLink	ConnexxAir
Programming adapter/cable	-

(1) Requires bilateral fitting

Fresh colors with warm metallic hues



Cosmic Blue & Rose Gold




Dark Granite & Silver



Snow White & Rose Gold

Styletto | Fitting Range



Styletto behind the ear

* Study conducted at the University of Northern Colorado, 2015, examined the effectiveness of the new features of primax by collecting and analyzing ongoing EEG data while subjects performed speech testing. For both primax features SpeechMaster and EchoShield, the objective brain behavior measures revealed a significant reduction in listening effort when the feature was activated.
2017 "First Fit Study" conducted at Hörzentrum Oldenburg examining the effect of own voice processing on spontaneous acceptance after first fit of hearing aids. Further details: www.signia-pro.com/first-fit-study.

** Two clinical studies have shown that binax provides better than normal hearing in certain demanding environments (University of Northern Colorado, 2014; Oldenburg Hörzentrum, 2013): Speech Reception Thresholds (SRT) in cocktail-party situations improved up to 2.9 dB for wearers with mild to moderate hearing loss using Carat binax or Pure binax hearing aids with narrow directionality, compared to people with normal hearing.

† Achieved IP68 rating per IEC 60529 standard.